

URETHRAL CATHETERISATION DEVICE (UCD®)

Video Demonstration & Feedback


We love to hear from you and support product adoption at your institution! Please scan the QR Code for product video demonstration and link to feedback form; or visit <https://scnv.io/dVOC>



Device Description

The innovative Urethrotech UCD® (Urethral Catheterisation Device) is composed of a 3-way Silicone Foley Catheter with integrated introducer guidewire for safe bladder access to prevent Catheterisation-associated Urethral Injury (CAUI). The integrated removable guidewire is made of Nitinol with hydrophilic coating to facilitate the atraumatic introduction of the catheter through the urethra into the bladder to negotiate any prostatic curve or tight external urethral sphincter. Also suitable for irrigation of the bladder after guidewire removal.

- >> Completely 100% Silicone 3-way Foley catheter
- >> Clear for easy visualisation of urine flow or debris
- >> Radiopaque blue tip and stripe for easy visualisation during scans
- >> Two large opposing eyes for efficient urine drainage
- >> Soft atraumatic integrated hydrophilic Nitinol guidewire
- >> Guidewire exits at tip of catheter for ideal catheter guidance
- >> Secure guidewire Luer-lock for lubrication of guidewire
- >> Attached plug ready to seal guidewire channel after guidewire removal
- >> Symmetrical balloon shape for added reliability
- >> Latex Free



Urethral Catheterisation Device (UCD®)

3-Way 100% Silicone Foley Catheter with integrated guidewire

PACKAGING						1box/each
Item No.	Size	Balloon	Length	NPC code	NHS-Codes Private CCSD-	Pkg.
UCD-0035-090-0016	16 Fr	10 mL	90cm	FST951	M478; M371 M4820	1/bx

Specifications

HPIS CODE	730_10_40_20
BALLOON SIZE (CC)	10
CATHETER MATERIAL	Silicone, 100% Silicone
COATING	No
FRENCH SIZE (FR)	16
LATEX FREE	Yes
PRODUCT STYLE	3-Way Foley Catheter, UCD
INTEGRATED GUIDEWIRE (CM)	90
SILVER	No
TEMPERATURE SENSING	No
UNSPSC	42142702

Unit of Measure	Conversion	Net Weight (g)	Volume (cubic m ³)	Shipping Dimensions L x W x H (cm)	GTIN-13
Box (BX)	1.0 Each	145	0.014742	52.5 x 10.8 x 2.6	5060639870004
Each (EA)	1.0 Each	145	0.014742	52.5 x 10.8 x 2.6	5060639870004

Product Information www.urethrotech.com



Contact

Product Training: training@urethrotech.com

Customer Service: sales@urethrotech.com

Direct Order: sales@urethrotech.com

Urethrotech Ltd
24 Garth Road, Kingston Upon Thames, Surrey KT2 5NY, UK
Tel: 0800 051 7219 www.urethrotech.com



Order via NHS Supply Chain (Code FST951)

